



Yes, you can.®

**Price List and Order Form**

Effective July, 2010

For ease of ordering, contact Customer Service Toll Free at: 0508-468-222 Fax: 0508-807-788

www.invacare.co.nz

Group \_\_\_\_\_ Date \_\_\_\_\_

Therapist/Email \_\_\_\_\_ Therapist Phone \_\_\_\_\_

Client \_\_\_\_\_ Therapist Fax \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

**POWERED SEATING BASE**

- TDXSP TDXSP with Formula CG Tilt w/ MK6i™ ..... 12375.00
- CGTE TDXSP CG Elevate to combine with Tilt ..... 14570.00
- CGTR TDXSP CG Recline to combine with Tilt ..... 16765.00
- CGTRE TDXSP CG Elevate/Recline with Tilt ..... 18960.00

**SEAT SELECTION**

**BACKS**

- BPO Back Posts Only ..... No Charge
- CURV Curved ..... 395.00

**SEAT WIDTH RANGE**

- 1216W 12" to 16" Seat Width Range (Adjustable) .... No Charge
- 1620W 16" to 20" Seat Width Range (Adjustable) ..... Standard
- 2024W 20" to 24" Seat Width Range (Adjustable) .... No Charge

**SEAT WIDTH SETTING**

- W12 12" Seat Width ..... No Charge
- W13 13" Seat Width ..... No Charge
- W14 14" Seat Width ..... No Charge
- W15 15" Seat Width ..... No Charge
- W16 16" Seat Width ..... No Charge
- W17 17" Seat Width ..... No Charge
- W18 18" Seat Width ..... No Charge
- W19 19" Seat Width ..... No Charge
- W20 20" Seat Width ..... No Charge
- W21 21" Seat Width ..... No Charge
- W22 22" Seat Width ..... No Charge
- W23 23" Seat Width ..... No Charge <sup>3</sup>
- W24 24" Seat Width ..... No Charge <sup>3</sup>

**SEAT DEPTH RANGE**

- 1218D 12" to 18" Seat Depth Range (Adjustable) .... No Charge
- 1619D 16" to 19" Seat Depth Range (Adjustable) ..... Standard
- 1922D 19" to 22" Seat Depth Range (Adjustable) .... No Charge

**SEAT DEPTH SETTING**

- D1212" Seat Depth ..... No Charge
- D1313" Seat Depth ..... No Charge
- D1414" Seat Depth ..... No Charge
- D1515" Seat Depth ..... No Charge
- D1616" Seat Depth ..... No Charge
- D1717" Seat Depth ..... No Charge
- D1818" Seat Depth ..... No Charge
- D1919" Seat Depth ..... No Charge
- D2020" Seat Depth ..... No Charge
- D2121" Seat Depth ..... No Charge
- D2222" Seat Depth ..... No Charge

**Base Price Includes**

- TDXSP for Formula CG Tilt with MK6i
- 136 Kg Weight Capacity (181kg Heavy Duty Option)
- 14" X 3" Pneumatic Tires ..... Standard
- 4 Pole Motors ..... Standard
- Flip Back Height Adjustable Arms ..... Standard
- MK690ACC Controller w/ Actuator Control ..... Standard
- 70NHD 70° Swingaway Footrests ..... Standard
- SPJ+ Joystick ..... Standard
- 1515 Push Button Style Seat Positioning Strap ..... Standard

**SYSTEM TYPE**

- CNVACTS Conventional Actuator ..... Standard
- SMTACTS Smart Actuator ..... 715.00

**BASE**

**USER WEIGHT LIMIT**

- U300 User Weight 136 kg ..... Standard <sup>1</sup>
- U400 User Weight 181 kg (Not for Elevate/Recline) .1095.00 <sup>2</sup>

**WHEEL OPTIONS (PAIR)**

- 1430-3 14" x 3" Pneumatic Tires ..... Standard
- 1431-3 14" x 3" with Foam Filled Inserts ..... 275.00
- Wheelchair Tie Down Brackets ..... 395.00

**FRAME FINISH**

**CONTEMPORARY**

- 61P Electric Red ..... Standard
- 62P Electric Blue ..... Standard

**SEAT**

**SEAT TILT**

- OSA 0° - 55° Seat Tilt ..... Standard

**FOOTNOTES**

<sup>1</sup> Package includes 4 Pole Motors and 5.8 mph top speed with <136 kg user weight limit. Package includes Quiet Stability Lock.  
<sup>2</sup> Package includes Heavy-Duty 4-pole motors and 5 mph top speed with <181 kg user weight limit.

**FOOTNOTES**

<sup>3</sup> Available only with 181 kg. package.

**NOTE: All specifications and dimensions are approximate.**

# Invacare® Formula™ CG Tilt/Recline/Elevate for TDX® SP

## Accessories and Options

### BACK HEIGHT

- TSBH20 Back Height 20" .....Standard
- TSBH21 Back Height 21" ..... No Charge
- TSBH22 Back Height 22" ..... No Charge
- TSBH23 Back Height 23" ..... No Charge
- TSBH24 Back Height 24" ..... No Charge
- TSBH25 Back Height 25" ..... No Charge
- TSBH26 Back Height 26" ..... No Charge

### ADJUSTABLE ARM TYPE (PAIR)

- 19F Reclining Full Length Adj. Hgt. Right..... No Charge
- 29F Reclining Desk Length Adj. Hgt. Left .....Standard

### STEALTH COMFORT HEADRESTS

- CP175 14" Multi Axis Headrest .....250.00
- CP175HW Removable & Adj Hardware .....250.00
- CP275 10" Multi Axis Headrest .....250.00
- CP275HW Removable & Adj Hardware .....250.00

## ELECTRONICS

### EXPANDABLE CONTROLLER SYSTEMS

- SPJ+ SPJ+ Joystick .....Standard
- MPJM6 MPJ™+ Multiple Drive Joystick .....750.00
- ASLRDYM6 MK6 Display Only (No Driver Control) ..... 825.00 <sup>4</sup>

### JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retractable Joystick Mnt SPJ+ .....435.00
- QLAM6 Quad Link Retractable Joystick Mnt MPJ+ .....435.00
- LEFT Left Handed Mounted Joystick ..... No Charge

### POWERED SEATING ACCESSORIES

- SANODE Single Act. Oper. through Driver Ctrl .....Standard <sup>5</sup>
- S4WSB Multiple Actuator Interface Box - must be Selected if more than one actuator function i.e. Tilt and Elevate Seat .....715.00
- FWT 4 Way Toggle ..... No Charge
- QPB 4 Push Buttons..... No Charge

## FRONT RIGGING OPTIONS

### FRONT RIGGING <sup>6</sup>

#### CENTER MOUNT <sup>7</sup>

- CTMT-F Center Mount Adjustable Knee Angle .....242.00 <sup>8</sup>
- CTMT-FJ Jr. Center Mount Adjustable Knee to Heel .....242.00 <sup>9</sup>

#### CENTER PIVOT STYLE (PAIR) <sup>10</sup>

- 70NHD 70° Swingaway Footrests .....Standard <sup>11</sup>
- 70STAPER 70° Tapered Footrests .....127.00 <sup>12</sup>
- AT5544 Manual Elevating Legrests.....352.00 <sup>13</sup>

### FOOTNOTES

- <sup>4</sup> Display Only, No Driver control.
- <sup>5</sup> Required for Tilt Operation through Driver Control – Single Actuator Systems Only. Not compatible with SPJ+ joysticks. Does not come standard with egg switch.
- <sup>6</sup> All front riggings include a standard footplate .
- <sup>7</sup> Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.
- <sup>8</sup> Extension range 12.5"-15".
- <sup>9</sup> Extension range 9.25"-12".
- <sup>10</sup> Provides up to 1 ½" of depth adjustment in ½" increments.
- <sup>11</sup> With seat widths of 21" to 24" and depth of 16", must pick footplate 1350. The AT5543 is not available..
- <sup>12</sup> Standard with impact guards. Not available with widths 19" to 24" and 16" deep seat.

### POWER ELEVATING LEGRESTS-ARTICULATING – INCLUDES ELECTRONICS

- MPCTMT Power Center Mount Foot Platform ..... 2195.00 <sup>14</sup>
- IFP Individual Footplates ..... 275.00 <sup>15</sup>
- MPCTMT-L Pwr Center Mnt Foot Platform Long..... 2195.00 <sup>14</sup>
- ELRPW Power ELR w/Articulation 20" Ext ..... 2195.00 <sup>16</sup>
- ELRPW22 Power ELR w/Articulation 22" Ext ..... 2195.00 <sup>17</sup>

### FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip-Up Footplates..... Standard <sup>18</sup>
- AT5543 Adj Angle Flip-Up Footplates ..... 205.00 <sup>19</sup>
- CTMTFB Footboard for Center Mount..... 165.00

### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap ..... 39.00
- 1600BK Heel Loops w/Ankle Straps) ..... 72.00
- IMPACTG Impact Guards ..... 72.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... 99.00

### MISCELLANEOUS ACCESSORIES

- 2ATP Perspex Tray – Swing away ..... 545.00

Motion Power Center Mount Knee to Heel			
STF Pick	STF at 0 deg	KTH w/IFP on MPCTMT	KTH with Footboard
Low	17.5"	10.5" - 14.5"	12" - 13.5" MPCTMT only
Med	18.5"	11.5" - 15.5"	12" - 14.5" MPCTMT only
High	19.5"	11.5" - 16.5"	12" - 15.5" MPCTMT or MPCTMT-L

### FOOTNOTES

- <sup>13</sup> Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.
- <sup>14</sup> Includes a standard, single footboard.
- <sup>15</sup> Individual footplate option for Power Center Mount. This option is only available on MPCTMT.
- <sup>16</sup> Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.
- <sup>17</sup> Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 15"-22". Sold as pair.
- <sup>18</sup> Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- <sup>19</sup> One piece footboard for Center Mount (CTMT-F) Riggings.