

# Invacare Action3 Junior - Mnl Chair - Category 1



Yes, you can.®

Price List and Order Form - Effective Feb, 2009

Invacare New Zealand , Ph: 0508 468 222, Fax: 0508 807 788, www.invacare.co.nz

Company Name:	Date:
Therapist:	Therapist Phone:
Therapist Email Address:	
Client:	Therapist fax:
Ship To Name and Address:	

**Includes as standard:** Push to lock brakes, positioning buckle belt, rear anti tippers, quick release axle on self propulsion, Desk flip up armrests

**x Please make your selection with a cross in the box**

<b>ACTION<sup>3</sup> JUNIOR EVOLUTIVE</b>			<b>BACKREST UPHOLSTERY TYPE</b>		
<input type="checkbox"/>	SELF PROPULSION	\$1,925.00	<input type="checkbox"/>	WIDTH & HEIGHT ADJUSTABLE NYLON SLING TYPE	standard
<input type="checkbox"/>	TRANSIT PROPULSION	\$1,925.00	<input type="checkbox"/>	TENSION ADJUSTABLE NYLON SLING TYPE	\$182.00
<b>SEAT</b>			<b>HANGERS/FOOTPLATES</b>		
<b>ADJUSTABLE SEAT WIDTH/DEPTH FRAME RANGE</b>			<b>HANGERS TYPE</b>		
<input type="checkbox"/>	8" - 11" Width Range x 8" - 10" Depth Range	n/c	<input type="checkbox"/>	90° STRAIGHT REMOVABLE HANGER; 12,5 TO 26 CM	standard
<input type="checkbox"/>	8" - 11" Width Range x 10" - 12" Depth Range	n/c	<input type="checkbox"/>	80° SHORT HANGER; 22 TO 32 CM	n/c
<input type="checkbox"/>	8" - 11" Width Range x 12" - 14" Depth Range	n/c	<input type="checkbox"/>	80° LONG HANGER; 30 TO 45 CM	n/c
<input type="checkbox"/>	12" - 15" Width Range x 11" - 13" Depth Range	n/c	<input type="checkbox"/>	ELEVATING LEGRESTS (0° TO 80°) INCLUDING CALFPAD); 27,5 TO 37,5 CM	\$352.00
<input type="checkbox"/>	12" - 15" Width Range x 13" - 15" Depth Range	n/c	<b>FOOTPLATES TYPE</b>		
<b>SEAT WIDTH SETTING</b>			<input type="checkbox"/>	FLIP UP FOOTPLATE	Standard
<input type="checkbox"/>	8" Seat Width	n/c	<input type="checkbox"/>	FLIP UP, ANGLE & DEPTH ADJUSTABLE FOOTPLATE	\$205.00
<input type="checkbox"/>	9" Seat Width	n/c	<b>ARMRESTS</b>		
<input type="checkbox"/>	10" Seat Width	n/c	<input type="checkbox"/>	REMOVABLE ADJ. 15 TO 25 CM (DELIVERED AT 18 CM) incl. Skirt Guard	Standard
<input type="checkbox"/>	11" Seat Width	n/c	<b>COLOURS</b>		
<input type="checkbox"/>	12" Seat Width	n/c	<input type="checkbox"/>	DARK RED	Std
<input type="checkbox"/>	13" Seat Width	n/c	<input type="checkbox"/>	SAPPHIRE BLUE	Std
<input type="checkbox"/>	14" Seat Width	n/c	<b>FRONT WHEELS</b>		
<input type="checkbox"/>	15" Seat Width	n/c	<input type="checkbox"/>	5" x 1"	Standard
<b>SEAT DEPTH SETTING</b>			<b>REAR WHEELS</b>		
<input type="checkbox"/>	8" Seat Depth	n/c	<input type="checkbox"/>	22" x 1.3/8" Pneumatic with Anodized Handrim	Standard
<input type="checkbox"/>	9" Seat Depth	n/c	<input type="checkbox"/>	20" x 1.3/8" Pneumatic with Anodized Handrim	\$149.00
<input type="checkbox"/>	10" Seat Depth	n/c	<input type="checkbox"/>	12" x 2" Pneumatic	Standard
<input type="checkbox"/>	11" Seat Depth	n/c	<input type="checkbox"/>	12" x 2" Solid	\$72.00
<input type="checkbox"/>	12" Seat Depth	n/c	<b>BRAKES TYPE</b>		
<input type="checkbox"/>	13" Seat Depth	n/c	<input type="checkbox"/>	PUSH TO LOCK BRAKES	Standard
<input type="checkbox"/>	14" Seat Depth	n/c	<b>SAFETY OPTIONS TYPE</b>		
<input type="checkbox"/>	15" Seat Depth	n/c	<input type="checkbox"/>	BUCKLE BELT	Standard
<b>SEAT TO FLOOR HEIGHT</b>			<input type="checkbox"/>	REAR ANTI-TIPPERS PAIR WITH TIPPING AID	Standard
<input type="checkbox"/>	16" Seat Height	n/c	<input type="checkbox"/>	BRAKE EXTENSION PAIR	\$39.00
<input type="checkbox"/>	17" Seat Height	n/c	<b>ACCESSORIES TYPE</b>		
<input type="checkbox"/>	18" Seat Height	n/c	<input type="checkbox"/>	SPOKE GUARD SMILEY	\$195.00
<b>SEAT UPHOLSTERY TYPE</b>			<input type="checkbox"/>	SPOKE GUARD SPIDER	\$195.00
<input type="checkbox"/>	NYLON WITH VELCRO STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE	Standard	<input type="checkbox"/>	SPOKE GUARD FLOWER	\$195.00
<b>BACKREST TYPE</b>			<input type="checkbox"/>	ADJUSTABLE HEELS STRAP	\$39.00
<input type="checkbox"/>	FIXED 0° BACKREST 3 HEIGHTS: 15"/16"/17" CM (DELIVERED AT 16")	n/c	<input type="checkbox"/>	ADJUSTABLE CALF STRAP	\$39.00
<input type="checkbox"/>	FIXED 7° BACKREST 3 HEIGHTS: 15"/16"/17" CM (DELIVERED AT 16")	Standard	<input type="checkbox"/>	ACRYLIC TRAY	\$250.00
<input type="checkbox"/>	FIXED 7° BACKREST 16" HEIGHT ADJUSTABLE PUSHHANDLE (6")	\$149.00	<input type="checkbox"/>	SWING AWAY ACRYLIC TRAY	\$435.00

\* All specs and dimensions are approximate. Non-standard items may incur a delay in delivery.