

- |                          |                        |
|--------------------------|------------------------|
| <input type="checkbox"/> | Purchase               |
| <input type="checkbox"/> | Trial                  |
| <input type="checkbox"/> | Quotation for purchase |
| <input type="checkbox"/> | Quotation for trial    |

Company Name: \_\_\_\_\_ Date: \_\_\_\_\_

Therapist: \_\_\_\_\_ Therapist Phone: \_\_\_\_\_

Client: \_\_\_\_\_ Therapist Fax: \_\_\_\_\_

Ship to Name and Address: \_\_\_\_\_

- Please make your selection with a cross in the box  
 Standard feature  
 Free of charge

CODE	DESCRIPTION		NZ Sell Price
DDM0012	küschall Compact	<input type="checkbox"/>	\$3,087.08
DDM0000	Special Construction ! Full description by fax or E-Mail !	<input type="checkbox"/>	

**SEATING CHOICE**

**FRONT FRAME**

DDM0211	Standard	<input checked="" type="checkbox"/>	
DDM0223	Hemiplegia frame, straight, left	<input type="checkbox"/>	\$110.13
DDM0210	Hemiplegia frame, straight, right	<input type="checkbox"/>	\$110.13
DDM0224	Hemiplegia frame, abducted, left	<input type="checkbox"/>	\$169.62
DDM0205	Hemiplegia frame, abduc., right	<input type="checkbox"/>	\$169.62

**SEAT WIDTH (SB) ! Double cross brace from SB44 cm onward !**

DDM0122	Seat width 32 cm	<input type="checkbox"/>	\$123.34
DDM0124	Seat width 34 cm	<input type="checkbox"/>	\$123.57
DDM0126	Seat width 36 cm	<input type="checkbox"/>	\$123.57
DDM0128	Seat width 38 cm	<input type="checkbox"/>	
DDM0130	Seat width 40 cm	<input type="checkbox"/>	
DDM0132	Seat width 42 cm	<input type="checkbox"/>	
DDM0134	Seat width 44 cm	<input type="checkbox"/>	\$19.09
DDM0136	Seat width 46 cm	<input type="checkbox"/>	\$21.65
DDM0138	Seat width 48 cm	<input type="checkbox"/>	\$15.04
DDM0140	Seat width 50 cm	<input type="checkbox"/>	\$143.77

**SEAT DEPTH (ST) ! Frame is adjusted to seat depth for maximum compactness !**

DDM0408	Seat depth 32 cm	<input type="checkbox"/>	
DDM0410	Seat depth 34 cm	<input type="checkbox"/>	
DDM0412	Seat depth 36 cm	<input type="checkbox"/>	
DDM0415	Seat depth 38 cm	<input type="checkbox"/>	
DDM0417	Seat depth 40 cm	<input checked="" type="checkbox"/>	
DDM0420	Seat depth 42 cm	<input type="checkbox"/>	
DDM0422	Seat depth 44 cm	<input type="checkbox"/>	
DDM0424	Seat depth 46 cm	<input type="checkbox"/>	\$144.30
DDM0427	Seat depth 48 cm	<input type="checkbox"/>	\$144.30
DDM0429	Seat depth 50 cm	<input type="checkbox"/>	\$144.30

**BACK HEIGHT (RH) !Foldable back not with RH 30-33 cm - Foldable back + angle adjustable back not with RH 30-40.5 cm!**

**Short Telescope DDM0597**


DDM0550	RH 30 cm	<input type="checkbox"/>	No Charge
DDM0503	RH 31.5 cm	<input type="checkbox"/>	No Charge
DDM0506	RH 33 cm	<input type="checkbox"/>	No Charge
DDM0509	RH 34.5 cm	<input type="checkbox"/>	No Charge
DDM0512	RH 36 cm	<input type="checkbox"/>	No Charge
DDM0514	RH 37.5 cm	<input type="checkbox"/>	No Charge
DDM0516	RH 39 cm	<input type="checkbox"/>	No Charge
DDM0518	RH 40.5 cm	<input type="checkbox"/>	No Charge
DDM0520	RH 42 cm	<input type="checkbox"/>	No Charge

**Medium telescope DDM0598**

DDM0516	RH 39 cm	<input type="checkbox"/>	\$1.34
DDM0518	RH 40.5 cm	<input type="checkbox"/>	\$1.34
DDM0520	RH 42 cm	<input checked="" type="checkbox"/>	\$1.34
DDM0522	RH 43.5 cm	<input type="checkbox"/>	\$1.34
DDM0524	RH 45 cm	<input type="checkbox"/>	\$1.34
DDM0526	RH 46.5 cm	<input type="checkbox"/>	\$1.34
DDM0528	RH 48 cm	<input type="checkbox"/>	No Charge
DDM0530	RH 49.5 cm	<input type="checkbox"/>	No Charge
DDM0532	RH 51 cm	<input type="checkbox"/>	No Charge

**Long Telescope DDM0599**

DDM0534	RH 52.5 cm	<input type="checkbox"/>	No Charge
DDM0536	RH 54 cm	<input type="checkbox"/>	No Charge
DDM0537	RH 55.5 cm	<input type="checkbox"/>	No Charge



CODE	DESCRIPTION		NZ sell price
<b>FRONT SEAT-TO-FLOOR HEIGHT (SHv)</b>			
<b>Standard front frame</b>			
DDM0844	SHv 42 cm	<input type="checkbox"/>	No Charge
DDM0845	SHv 43 cm	<input type="checkbox"/>	No Charge
DDM0846	SHv 44 cm	<input type="checkbox"/>	No Charge
DDM0847	SHv 45 cm	<input type="checkbox"/>	No Charge
DDM0848	SHv 46 cm	<input type="checkbox"/>	No Charge
DDM0849	SHv 47 cm	<input type="checkbox"/>	No Charge
DDM0850	SHv 48 cm	<input type="checkbox"/>	No Charge
DDM0851	SHv 49 cm	<input type="checkbox"/>	\$3.67
DDM0852	SHv 50 cm	<input checked="" type="checkbox"/>	\$3.67
DDM0853	SHv 51 cm	<input type="checkbox"/>	\$3.29
DDM0854	SHv 52 cm	<input type="checkbox"/>	\$32.38
DDM0855	SHv 53 cm	<input type="checkbox"/>	
DDM0856	SHv 54 cm	<input type="checkbox"/>	
DDM0857	SHv 55 cm	<input type="checkbox"/>	
<b>Hemiplegia front frame</b>			
DDM0839	SHv 37 cm	<input type="checkbox"/>	
DDM0840	SHv 38 cm	<input type="checkbox"/>	
DDM0841	SHv 39 cm	<input type="checkbox"/>	
DDM0842	SHv 40 cm	<input type="checkbox"/>	
DDM0843	SHv 41 cm	<input type="checkbox"/>	
DDM0844	SHv 42 cm	<input checked="" type="checkbox"/>	
DDM0845	SHv 43 cm	<input type="checkbox"/>	
DDM0846	SHv 44 cm	<input type="checkbox"/>	
DDM0847	SHv 45 cm	<input type="checkbox"/>	
DDM0848	SHv 46 cm	<input type="checkbox"/>	
DDM0849	SHv 47 cm	<input type="checkbox"/>	
DDM0850	SHv 48 cm	<input type="checkbox"/>	
DDM0851	SHv 49 cm	<input type="checkbox"/>	
DDM0852	SHv 50 cm	<input type="checkbox"/>	
<b>REAR SEAT-TO-FLOOR HEIGHT(SHh) ! SHh 36 to 38 cm only for hemiplegic option !</b>			
DDM0941	SHh 36 cm	<input type="checkbox"/>	
DDM0942	SHh 37 cm	<input type="checkbox"/>	
DDM0943	SHh 38 cm	<input type="checkbox"/>	
DDM0944	SHh 39 cm	<input type="checkbox"/>	\$3.54
DDM0945	SHh 40 cm	<input type="checkbox"/>	\$3.54
DDM0946	SHh 41 cm	<input type="checkbox"/>	\$3.54
DDM0947	SHh 42 cm	<input type="checkbox"/>	\$3.54
DDM0948	SHh 43 cm	<input type="checkbox"/>	\$3.54
DDM0949	SHh 44 cm	<input type="checkbox"/>	
DDM0950	SHh 45 cm	<input checked="" type="checkbox"/>	
DDM0951	SHh 46 cm	<input type="checkbox"/>	
DDM0952	SHh 47 cm	<input type="checkbox"/>	
DDM0953	SHh 48 cm	<input type="checkbox"/>	
DDM0954	SHh 49 cm	<input type="checkbox"/>	No Charge
DDM0955	SHh 50 cm	<input type="checkbox"/>	No Charge
<b>KNEE-TO-HEEL LENGTH ! UL ≤ (SHv - 2cm) ! ! UL51 and 52 cm not with legrest 80° !</b>			
<b>UL range for footplate high-mounted DDM1589</b>			
DDM1661	UL 20 cm	<input type="checkbox"/>	
DDM1662	UL 21 cm	<input type="checkbox"/>	
DDM1663	UL 22 cm	<input type="checkbox"/>	
DDM1664	UL 23 cm	<input type="checkbox"/>	
DDM1665	UL 24 cm	<input type="checkbox"/>	
DDM1666	UL 25 cm	<input type="checkbox"/>	
DDM1667	UL 26 cm	<input type="checkbox"/>	
DDM1668	UL 27 cm	<input type="checkbox"/>	
DDM1669	UL 28 cm	<input type="checkbox"/>	
DDM1670	UL 29 cm	<input type="checkbox"/>	
DDM1671	UL 30 cm	<input type="checkbox"/>	
DDM1672	UL 31 cm	<input type="checkbox"/>	
<b>Short UL telescope DDM1590</b>			
DDM1673	UL 32 cm	<input type="checkbox"/>	No Charge
DDM1674	UL 33 cm	<input type="checkbox"/>	No Charge
DDM1675	UL 34 cm	<input type="checkbox"/>	No Charge
DDM1676	UL 35 cm	<input type="checkbox"/>	No Charge
DDM1677	UL 36 cm	<input type="checkbox"/>	No Charge
DDM1678	UL 37 cm	<input type="checkbox"/>	No Charge
DDM1679	UL 38 cm	<input type="checkbox"/>	No Charge
DDM1680	UL 39 cm	<input type="checkbox"/>	No Charge
DDM1681	UL 40 cm	<input type="checkbox"/>	No Charge
<b>Medium UL telescope DDM1591</b>			
DDM1679	UL 38 cm	<input type="checkbox"/>	
DDM1680	UL 39 cm	<input type="checkbox"/>	
DDM1681	UL 40 cm	<input type="checkbox"/>	
DDM1682	UL 41 cm	<input type="checkbox"/>	
DDM1683	UL 42 cm	<input type="checkbox"/>	
DDM1684	UL 43 cm	<input checked="" type="checkbox"/>	
DDM1685	UL 44 cm	<input type="checkbox"/>	
DDM1686	UL 45 cm	<input type="checkbox"/>	
DDM1687	UL 46 cm	<input type="checkbox"/>	

**SHv / Castor Size**

\* Compact Standard

3" : 42-51 cm  
 4" : 44-52 cm  
 5" : 47-53 cm  
 6" : 49-54 cm  
 7" : 52-55 cm

**SHv / Castor Size**


\* Compact Hemiplegic

3" : 37-45 cm  
 4" : 39-46 cm  
 5" : 41-47 cm  
 6" : 44-49 cm  
 7" : 46-50 cm











**SHh / Rear wheel size**







22" : 37-45 cm  
 24" : 39-48 cm  
 25" : 41-50 cm  
 26" : 42-50 cm










Customer UL must be measured taking away the cushion thickness








**UL = ULpat - SK**

CODE	DESCRIPTION		NZ sell price	
<b>Long UL telescope DDM1592</b>				
DDM1685	UL 44 cm	<input type="checkbox"/>	\$0.68	
DDM1686	UL 45 cm	<input type="checkbox"/>	\$0.68	
DDM1687	UL 46 cm	<input type="checkbox"/>	\$0.68	
DDM1688	UL 47 cm	<input type="checkbox"/>	\$0.68	
DDM1689	UL 48 cm	<input type="checkbox"/>	\$0.68	
DDM1690	UL 49 cm	<input type="checkbox"/>		
DDM1691	UL 50 cm	<input type="checkbox"/>		
DDM1692	UL 51 cm	<input type="checkbox"/>		
DDM1693	UL 52 cm	<input type="checkbox"/>		
<b>BACKREST CHOICE</b>				
<b>BACK UPHOLSTERY</b>				
DDM0700	Standard, light upholstery	<input type="checkbox"/>		
DDM0701	Tension velcro adjustable back	<input type="checkbox"/>	\$20.58	
DDM0735	Aero back, tension velcro adjustable	<input type="checkbox"/>	\$245.97	
<b>PUSH HANDLES ! if no option selected, standard push handles will be assembled !</b>				
DDM0703	Height adjustable rear set ! not with standard, back upholstery !	<input type="checkbox"/>	\$125.24	
DDM0733	Height adjustable, integrated	<input type="checkbox"/>	\$180.35	
DDM0704	Foldable push handles !not with drum brakes, not with height adj. integr. push handles	<input type="checkbox"/>	\$102.91	
<b>BACK TUBE</b>				
DDM0705	Lumbar back	<input type="checkbox"/>	\$5.37	
DDM0734	Foldable back !not with height adj. Integrated push handles and std back upholstery !	<input type="checkbox"/>	\$217.24	
<b>BACK ANGLE ADJUSTABLE</b>				
DDM0706	Back Angle -12°	<input type="checkbox"/>	\$162.08	
DDM0707	Back Angle -8°	<input type="checkbox"/>	\$162.08	
DDM0718	Back Angle -4°	<input type="checkbox"/>	\$162.08	
DDM0709	Back Angle 0°	<input type="checkbox"/>	\$162.08	
DDM0719	Back Angle +4°	<input type="checkbox"/>	\$162.08	
DDM0711	Back Angle +8°	<input type="checkbox"/>	\$162.08	
DDM0712	Back Angle +12°	<input type="checkbox"/>	\$162.08	
<b>SEATING</b>				
<b>SEAT UPHOLSTERY</b>				
DDM0611	Standard adjustable seat upholstery	<input type="checkbox"/>		
DDM0613	Seat upholstery with 1 central Velcro band (soft band)	<input type="checkbox"/>	\$1.32	
DDM0614	Seat upholstery with 2 side Velcro bands rough band	<input type="checkbox"/>	\$1.29	
<b>FOOTREST</b>				
<b>LEGREST ! 80° legrest not with 6" and 7" front wheels</b>				
DDM1540	Legrest, angle open, 70°	<input type="checkbox"/>		
DDM1541	Legrest, angle closed, 80°	<input type="checkbox"/>		
DDM1528	Legrest height adjustable left	<input type="checkbox"/>	\$827.85	
DDM1529	Legrest height adjustable right	<input type="checkbox"/>	\$827.85	
DDM1538	Legrest height adjustable both sides ! Not possible with SB 36 cm	<input type="checkbox"/>	\$1,655.70	
<b>FOOTREST ! not with a height adjustable legrest both sides - Rear set not with 7" castor !</b>				
<b>! 80° legrest not with 5" front wheels in combination with footplates set to back and seat angle smaller than 5 cm</b>				
DDM1562	2 pc plastic set to front	<input type="checkbox"/>		
DDM1563	2 pc plastic set to back	<input type="checkbox"/>		
DDM1564	2 piece, adj set to front	<input type="checkbox"/>	\$136.33	
DDM1565	2 piece, adj set to back	<input type="checkbox"/>	\$136.33	
DDM1566	1 piece, adj to front	<input type="checkbox"/>	\$243.70	
DDM1567	1 piece, adj set to back	<input type="checkbox"/>	\$243.70	
DDM1556	Footrest upmounted	<input type="checkbox"/>	\$342.94	
<b>FOOT AND LEG SUPPORT ! Loops only with 2 piece footplate or high mounted !</b>				
DDM1534	Calf strap	<input type="checkbox"/>		
DDM1530	Heel loops	<input type="checkbox"/>	\$18.76	
DDM1536	Instep loops	<input type="checkbox"/>	\$13.14	
<b>SIDE GUARDS</b>				
<b>SIDE GUARDS ! not with siderests !</b>				
DDM1400	Clothe guards	<input type="checkbox"/>		
DDM1402	Mud guards, fixed ! Not possible with rear wheel extension	<input type="checkbox"/>	\$65.70	
DDM1403	Mud guards removable ! Ads 2 cm to the total width of the wheelchair	<input type="checkbox"/>	\$269.80	
DDM1405	Mud guards carbon, fixed ! Not possible with rear wheel extension	<input type="checkbox"/>	\$154.05	
DDM1408	Mud guards carbon, removable ! Ads 2 cm to the total width of the wheelchair	<input type="checkbox"/>	\$293.42	
<b>SIDEREST ! included clothe guards ! Ads 2 cm to the total width of the wheelchair</b>				
DDM1417	Siderest flip-back removable, short	<input type="checkbox"/>	\$348.00	
DDM1420	Siderest flip-back removable, long	<input type="checkbox"/>	\$349.01	
DDM1416	Siderest height adjustable (with tools), short	<input type="checkbox"/>	\$245.16	
DDM1421	Siderest height adjustable (with tools), long	<input type="checkbox"/>	\$237.09	
DDM1430	Siderest flip-back height adjustable removable, short	<input type="checkbox"/>	\$271.04	
DDM1431	Siderest flip-back height adjustable removable, long	<input type="checkbox"/>	\$271.04	
<b>ARMREST</b>				
DDM1411	Tubular armrest, height adjustable, removable, swing-away, comfort padded	<input type="checkbox"/>	\$135.14	

CODE	DESCRIPTION		NZ sell price	
<b>SEATING OPTIONS</b>				
<b>SEAT CUSHION ! Vicair available for SB38, 40, 42, 44 !</b>				
DDM0600	Standard seat cushion	<input type="checkbox"/>	\$29.77	  0600                      0602
DDM0650	Vicair Comfair	<input type="checkbox"/>	\$605.90	
DDM0651	Vicair Incontinent	<input type="checkbox"/>	\$617.62	
DDM0601	Extra Standard cushion cover	<input type="checkbox"/>	\$25.11	
DDM0655	Extra Comfair cushion cover	<input type="checkbox"/>	\$76.99	
DDM0656	Extra Incontinent cushion cover	<input type="checkbox"/>	\$88.71	
DDM0602	Seat pouch for Standard cushion	<input type="checkbox"/>	\$46.18	
<b>FRAME COLOUR ! A colour difference of 7% might be possible between the front frame, back frame and cross-brace (Official ISO norm) !</b>				
DDM0300	Glossy-Black	<input type="checkbox"/>		
DDM0301	Bronze-Black	<input type="checkbox"/>	\$8.58	
DDM0302	Royal Blue	<input type="checkbox"/>	\$22.28	
DDM0303	Royal Red	<input type="checkbox"/>	\$22.28	
DDM0350	Mystic Purple	<input type="checkbox"/>	\$44.56	
DDM0315	Silver Grey	<input type="checkbox"/>	\$39.62	
DDM0351	Xanthic Yellow	<input type="checkbox"/>	\$6.23	
DDM0352	Cyber Green	<input type="checkbox"/>	\$43.70	
<b>COLOUR-POWER-SET ! small shade differences are possible ! ! rear wheel extension not in CPS !</b>				
DDM0330	Colour-Power-Set techno-anthracite	<input type="checkbox"/>	\$480.61	<b>CPS-Set :</b> Starec-hub, forks, adapter plates, rear wheels
DDM0331	Colour-Power-Set lemon	<input type="checkbox"/>	\$438.89	
DDM0332	Colour-Power-Set apricot	<input type="checkbox"/>	\$469.29	
DDM0333	Colour-Power-Set purple	<input type="checkbox"/>	\$469.29	
DDM0334	Colour-Power-Set b-blue	<input type="checkbox"/>	\$438.89	
<b>FRONT WHEEL ! 6" &amp; 7" not with legrest 80° !</b>				
DDM0800	3" Supersport	<input type="checkbox"/>	\$5.39	 0800    0801    0810   0811    0812    0820/23
DDM0801	4" Everyday	<input type="checkbox"/>		
DDM0802	4" Starec colourless	<input type="checkbox"/>	\$60.94	
DDM0803	4" Starec techno-anthracite	<input type="checkbox"/>	\$60.35	
DDM0804	4" Starec lemon	<input type="checkbox"/>	\$63.34	
DDM0805	4" Starec apricot	<input type="checkbox"/>	\$63.49	
DDM0806	4" Starec violett	<input type="checkbox"/>	\$63.49	
DDM0807	4" Starec b-blue	<input type="checkbox"/>	\$63.49	
DDM0809	4" Starec black	<input type="checkbox"/>	\$62.66	
DDM0824	4" Skater Wheel	<input type="checkbox"/>	\$3.09	
DDM0810	5" Light Wheel	<input type="checkbox"/>		
DDM0811	5" Softroll	<input checked="" type="checkbox"/>	\$2.28	
DDM0812	5" Starec colourless	<input type="checkbox"/>	\$65.01	
DDM0813	5" Starec techno anthracite	<input type="checkbox"/>	\$65.87	
DDM0814	5" Starec lemon	<input type="checkbox"/>	\$67.37	
DDM0815	5" Starec apricot	<input type="checkbox"/>	\$67.37	
DDM0816	5" Starec violett	<input type="checkbox"/>	\$67.37	
DDM0817	5" Starec b-blue	<input type="checkbox"/>	\$67.37	
DDM0819	5" Starec black	<input type="checkbox"/>	\$65.54	
DDM0820	6" Light wheel	<input type="checkbox"/>	-\$3.82	
DDM0821	6" Pneumatic Wheel	<input type="checkbox"/>	\$14.71	
DDM0822	7" Light Wheel	<input type="checkbox"/>	\$28.71	
DDM0823	7" Pneumatic Wheel	<input type="checkbox"/>	\$19.32	
<b>FRONT WHEEL ATTACHEMENT</b>				
DDM0867	Standard (Flow System)	<input checked="" type="checkbox"/>		
DDM0868	Bearing blocks (bringing the front wheels 2 cm outwards)	<input type="checkbox"/>		
<b>FRONT WHEEL EXTENSION ! Only possible with the Flow system</b>				
DDM0869	Front wheel extension	<input type="checkbox"/>		 0869
<b>QUICK RELEASE FORK ! Not possible with 7" Pneumatic front wheels</b>				
DDM0870	Quick release	<input type="checkbox"/>	\$272.81	
<b>2nd FRONT WHEEL SET FOR QUICK RELEASE</b>				
DDM0871	3" Supersport	<input type="checkbox"/>	\$228.71	
DDM0872	4" Everyday	<input type="checkbox"/>	\$228.71	
DDM0873	4" Starec colourless	<input type="checkbox"/>	\$228.71	
DDM0874	4" Starec techno-anthracite	<input type="checkbox"/>	\$289.06	
DDM0875	4" Starec lemon	<input type="checkbox"/>	\$292.03	
DDM0876	4" Starec apricot	<input type="checkbox"/>	\$292.20	
DDM0877	4" Starec violett	<input type="checkbox"/>	\$292.20	
DDM0878	4" Starec b-blue	<input type="checkbox"/>	\$292.20	
DDM0880	4" Starec black	<input type="checkbox"/>	\$291.34	
DDM0881	5" Light Wheel	<input type="checkbox"/>	\$228.71	
DDM0882	5" Softroll	<input type="checkbox"/>	\$226.43	
DDM0883	5" Starec colourless	<input type="checkbox"/>	\$293.72	
DDM0884	5" Starec techno anthracite	<input type="checkbox"/>	\$296.86	
DDM0885	5" Starec lemon	<input type="checkbox"/>	\$296.08	
DDM0886	5" Starec apricot	<input type="checkbox"/>	\$296.08	
DDM0887	5" Starec violett	<input type="checkbox"/>	\$296.08	
DDM0888	5" Starec b-blue	<input type="checkbox"/>	\$296.08	
DDM0890	5" Starec black	<input type="checkbox"/>	\$294.23	
DDM0891	6" Light wheel	<input type="checkbox"/>	\$228.33	
DDM0892	6" Pneumatic	<input type="checkbox"/>	\$254.94	
DDM0893	7" Light Wheel	<input type="checkbox"/>	\$264.73	

CODE	DESCRIPTION		NZ sell price	
<b>REAR WHEEL CHOICE</b>				
<b>REAR WHEEL FIXATION</b>				
DDM0250	Standard fixation	◆		
DDM0249	Rear wheel extension	□	\$58.25	
<b>REAR WHEEL POSITION ! When rear wheel extension is selected, position 5 and 6 only possible</b>				
DDM0263	1st position ! not with drum brakes ! ! Not with scissor brakes !	□	\$3.57	
DDM0264	2nd position ! Not with scissor brakes !	□	\$3.57	Std axle
DDM0265	3rd position ! Not with scissor brakes !	◆		
DDM0266	4th position ! 4th position not with stoplock brakes	□		Positions
DDM0267	5th position Not possible with mud guard fix	□	\$58.25	Rear wheel extension
DDM0268	6th position ! Not possible with mud guard fix	□	\$58.25	
<b>QUAD QUICK RELEASE AXLE</b>				
DDM0909	Quad quick release axle ! Not possible with stoplock brake !	□	\$45.67	
<b>REAR WHEEL SIZE</b>				
DDM0901	22"	□	\$249.77	
DDM0902	24"	□		
DDM0925	25" ! Only with High Performance wheel !	□	\$460.71	
DDM0903	26"	□	\$404.30	
<b>REAR WHEEL CAMBER</b>				
DDM0282	Camber 1°	◆		
DDM0283	Camber 4° (not with 26" rear wheels or performance brakes)	□	\$34.53	
<b>REAR WHEEL TYPE ! Drum brakes are delivered with specific wheels</b>				
DDM0904	24" Spider-Wheel	□	\$247.77	
DDM0914	22", 24", 26" küschall Standard	◆		
DDM0915	22", 24", 26" küschall Starec (Hub in CPS if selected)	□	\$20.23	
DDM0908	24" High Performance (Starec hub)	□	\$329.77	
DDM0918	25" High Performance (Starec hub) ! Only with "Marathon Evolution" light wheel tyre and light alu handrim	□	\$460.71	
DDM0913	24" Mountain Bike (complete wheel)	□	\$365.14	
DDM0917	24" Spinergy wheel	□	\$551.65	
<b>TYPES OF SPOKES ! only for küschall Starec wheel !</b>				
DDM0910	Cross Spoked	□		
DDM0911	Radial Spoked	□	\$388.23	
<b>REAR WHEEL TYRE ! Not for MB wheel ! ! 24" High Performance wheel only possible with light wheel, profil tyre and Marathon Evolution tyre ! 25" High Performance only available with "Marathon Evolution" tyre</b>				
DDM1001	22", 24", 26" Light wheel 7 bar	◆	\$30.51	
DDM1002	22", 24" Profil tyre 7 bar	□	\$153.90	
DDM1003	22", 24", 26" Sports hall collé	□	\$36.43	
DDM1004	24" Solid tyre, profil, grey	□	\$34.84	
DDM1005	24" Solid tyre (KIK) grey	□	\$34.84	
DDM1006	24" Solid tyre (KIK), black	□	\$36.43	
DDM1007	22", 24" Solid, light wheel	□	\$55.65	
DDM1010	24", 25" "Marathon Evolution" tyre from Schwalbe (puncture protected thanks to kevlar layer)	□	\$55.65	
<b>ANTI-FLAT</b>				
DDM1020	Anti-flat tyre tube protection	□	\$79.77	
<b>HANDRIM ! With 24" High Perf. &amp; Spinergy : alu light is mounted titanium &amp; Supergripp possible ! ! With 25" High performance wheels, only alu light possible ! ! With MB wheels, alu. is mounted, Supergripp possible ! Spider not with steel, Maxgrepp, alu light &amp; Tetra</b>				
DDM1100	22", 24", 26" Alu, distance 3 cm	◆		
DDM1101	22", 24", 26" Alu, distance 4 cm	□	\$30.61	
DDM1102	22", 24", 26" Supergripp, dist. 3 cm	□	\$128.10	
DDM1103	22", 24", 26" Supergripp, dist. 4 cm	□	\$129.62	
DDM1104	24", 26" Steel, dist. 3 cm	□	\$74.05	
DDM1105	24", 26" Steel, dist. 4 cm	□	\$74.20	
DDM1106	24" Titanium, dist. 3 cm	□	\$281.72	
DDM1107	24" Titanium, dist. 4 cm	□	\$281.87	
DDM1109	24" Maxgrepp	□	\$266.00	
DDM1110	24", 26" Tetra 2000	□	\$265.97	
DDM1113	24", 25" Light Alu, dist. 3 cm	□	\$35.62	
DDM1114	24", 25" Light Alu, dist. 4 cm	□	\$35.62	
<b>SPOKE GUARDS</b>				
<b>SPOKE GUARDS ! Not with Spider, Spinergy &amp; High Performance - For M.Bike: only 22" spokes guards ! ! Not possible with drum brakes</b>				
DDM1200	22", 24", 26" Transparent	□	\$87.01	
DDM1202	küsch (24")	□	\$115.82	
DDM1203	Nimbo (24")	□	\$126.38	
DDM1204	Boa-Boa (24")	□	\$141.44	
DDM1206	Cosmos Brownie (22",24")	□	\$137.49	
DDM1207	"Original" Carbon (24")	□	\$336.53	

SAFETY			NZ sell price	
<b>BRAKES</b>				
DDM1300	Push brake, bent	◆		
DDM1301	Push brake lever straight	□	\$2.20	
DDM1302	Performance brake ! Not with 4° camber and not with drum brakes !	□	\$59.27	
DDM1303	Scissor brake !Only possible with 4th rear wheel position with ST 40 cm or more	□	\$6.61	
DDM1306	Pull brake	□	\$6.61	
<b>DRUM BRAKE</b> ! Not with Quad release, 4° camber & foldable push handles, spoke guards ! Specific rear wheel mounted with a black adapterplate				
DDM1304	Attendant drum brake	□	\$400.20	
<b>STOPLOCK BRAKE</b> ! Not with Quad Axle, Drum brake, Spider, High Perf., M.Bike, Starec, Spinergy, 4° camber , 4th position !				
DDM1305	Stoplock brake	□	\$498.20	
<b>BRAKE EXTENSION LEVER</b> ! Only for brake lever straight or Pull brake !				
DDM1310	Brake extension lever	□	\$65.27	
<b>ANTI TIPPER</b>				
DDM1702	Antitipper right	□	\$165.04	
DDM1703	Antitipper left	□	\$165.04	
DDM1704	Antitipper both sides	□	\$312.10	
<b>SEAT BELT</b>				
DDM1730	Seat belt	□	\$15.34	
<b>PASSIVE LIGHTS</b>				
DDM1710	Passive lights assembled (mounted on spokes)	□		
DDM1711	Passive lights enclosed	□		
<b>ACCESSORIES</b>				
<b>LATERAL SUPPORT</b> ! Not with standard back light upholstery ! Not with Siderest flip-back removable and not with Siderest flip-back removable height adjustable				
DDM1731	Lateral Support, both sides	□	\$774.68	
DDM1732	Lateral Support, right	□	\$397.72	
DDM1733	Lateral Support, left	□	\$397.72	
<b>THERAPY TABLE</b> ! Only with siderests height adjsutable, removable short or long				
DDM1760	Table adjustable	□	\$259.87	
<b>TIPPER AID</b> ! Not on same side as cane holder !				
DDM1708	Tipper-Aid, right	□	\$64.46	
DDM1709	Tipper-Aid, left	□	\$64.46	
DDM1718	Tipper-Aid, both sides	□	\$130.89	
<b>CANE HOLDER</b> ! Not on the same side as tipper aid !				
DDM1706	Cane holder, right	□	\$93.85	
DDM1707	Cane holder, left	□	\$93.85	
<b>TRANSIT WHEELS 3"</b> ! Max. 2 adaptations per side !				
DDM1700	Transit wheels 3"	□	\$164.41	
<b>MANUAL AIR PUMP</b>				
DDM1715	Air pump	□		
<b>AIRMAN</b>				
DDM1716	AirMan	□	\$85.42	
<b>TOOL KIT</b>				
DDM1717	Tool kit	□		
<b>WHEEL BAG</b>				
DDM1728	Wheel bag (for 24" & 26" rear wheel)	□	\$43.87	
<b>WHEEL BAG &amp; WHEELS KITS</b>				
DDM1734	24" Mountain bike wheels with bag	□	\$342.00	
DDM1735	24" High Performance wheels with bag (handrims in Aluminium)	□	\$261.27	
DDM1737	25" High Performance wheels with bag (handrims in Aluminium)	□	not available	
<b>FREIGHT &amp; PRICE</b>				
FREIGHT PRICE				
TOTAL PRICE				