

Invacare New Zealand
 4 Westfield Place,
 PO Box 62124
 Mt Wellington
 Ph: 0508 468 222
 Fax: 0508 807 788



Purchase
 Trial
 Quotation for purchase
 Quotation for trial

Yes, you can.®

Invacare M71 Non Tilt

Max user weight 136kg

ACCESSIBLE PRICE LIST AND ORDER FORM - EFFECTIVE 1 SEPTEMBER 2006

Company Name: _____ Date: _____

Therapist: _____ Therapist Phone: _____

Client: _____ Therapist fax _____

Ship To Name and Address: _____

Includes as standard: Angle adjustable flip up metal footplates, 70 deg swing away legrest hangers, Calf Strap, Seat pan, Buckle seat belt, Height adjustable desk armrests, Straight 45cm high back canes, 40amp batteries, 10 1/2" x 3 1/2" Pneumatic tyres, True centre wheel drive technology, Dynamic shark electronics.

Please make your selection with a cross in the box

X	Description	Price
SEAT WIDTH		
	16" Seat width	n/c
	18" Seat width	Standard
	20" Seat width	\$490.00
ADJUSTABLE SEAT WIDTH RANGE		
	16" – 20" Seat depth range	Standard
SEAT DEPTH		
	16" Seat depth	Standard
	18" Seat depth	n/c
	20" Seat depth	\$490.00
ADJUSTABLE SEAT DEPTH RANGE		
	16" – 19" Seat depth range	Standard
	19" – 20" Seat depth range	n/c
TILT BACK OPTION BACK TYPE		
	Posts only	Standard
	Curved back	\$480.00
	Straight back canes (angle adj from base)	Standard
	Sectional rear uprights 10 deg angle	\$120.00
SEAT BACK HEIGHT		
	40cm (16") Back height	n/c
	43cm (17") Back height	n/c
	45cm (18") Back height	Standard
	48cm (19") Back height	n/c
	50cm (20") Back height	n/c
ARM TYPE		
	Desk length adjustable height	Standard
	Full length adjustable height	\$195.00
FOOTBOARD/LEGREST (HEMI)		
	Swingaway footrests 13 1/2" – 19"	Standard
	Swingaway ELR 16"-9 3/4" (pair)	\$465.00
	Smartleg Art. Legrests 17 1/4"-21 1/2" (pair)	\$495.00
	Footboard	\$375.00

X	Description	Price
FOOTPLATE OPTIONS		
	Angle adjustable flip up metal footplates	Standard
	Aluminium Footplates	\$90.00
	Extra Large Aluminium Footplates	\$140.00
FOOTREST/LEGREST ACCESSORIES		
	Impact guards (pair)	\$90.00
	Longer pivot and slide tube 2" (pair)	\$95.00
	Longer pivot and slide tube 4" (pair)	\$95.00
MISCELLANEOUS ACCESSORIES		
	Swing away joystick	\$330.00
	Oxygen holder	\$215.00
	Transport tie down brackets	\$380.00
	Lap tray clear with cut	\$295.00
SHROUD COLOURS		
	Deep blue pearl	Standard
	Candy red pearl	n/c
NON PREFERRED SUPPLY SEATING OPTIONS		
CUSHIONS		
	Comfort-Mate Extra Cushion	
	Infinity Gentle FloGel	
	Infinity Gentle ViscoFoam	
	Infinity Max FloGel	
	Infinity Max ViscoFoam	
	Ulti-Mate Cushion	
BACKS		
	Curved Back	
	Infinity Dual Flex	
	Personal Back 10	
	Personal Back 10 Plus	