

# Invacare® Spree GT - Pediatric Growable Tilt-in-Space Wheelchair



**Price List and Order Form** Effective October 18, 2004 Revised July 14, 2006

For ease of ordering, contact Customer Service Toll Free at:  
800-333-6900 Fax: 800-678-4682  
To order Literature: 800-828-6282

www.invacare.com

**Yes, you can.®**

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone # \_\_\_\_\_ P.O.# \_\_\_\_\_ Date \_\_\_\_\_

Ship To Name & Address \_\_\_\_\_

### Special Note

**HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.**

## FRAME SELECTIONS

### FRAME TYPE

E1233  SPREEGT Spree GT Tilt-in-Space Wheelchair .....2,754.00

### BACK TYPE

E1226  RECLINER Reclining Back.....810.00  
E0966  RECEXT Reclining Headrest Pillow and Extension.....75.00

### SEAT WIDTH RANGE

SWR1216 12"-16" Seat Width Range.....No Charge  
 SWR1418 14"-18" Seat Width Range.....No Charge

### FRAME WIDTH

W12 12" Wide .....No Charge  
 W13 13" Wide .....No Charge  
 W14 14" Wide .....No Charge  
 W15 15" Wide .....No Charge  
 W16 16" Wide .....No Charge  
 W17 17" Wide .....No Charge  
 W18 18" Wide .....No Charge

### SEAT FRAME LENGTH

SDR1216 12"-16" Seat Depth Range.....No Charge  
 SDR1519 15"-19" Seat Depth Range.....No Charge

### SEAT DEPTH

D12 12" Seat Depth .....No Charge  
 D13 13" Seat Depth .....No Charge  
 D14 14" Seat Depth .....No Charge  
 D15 15" Seat Depth .....No Charge  
 D16 16" Seat Depth .....No Charge  
 D17 17" Seat Depth .....No Charge  
 D18 18" Seat Depth .....No Charge  
 D19 19" Seat Depth .....No Charge

### SEAT STYLE

FLSEAT Flush Mounted Aluminum Seat Pan .....95.00  
 DPSEAT Drop Mounted Aluminum Seat Pan .....95.00

## TRANSPORT/VENT TRAY SELECTIONS

### TRANSPORT TIE DOWN

<sup>3,4</sup>  TRRO2 Transport Ready Option-Personal Back/Cushion .....No Charge  
 TRBKTS2 Transport Brackets .....No Charge  
 TROMIT Less Transport Ready Option.....Less 150.00

### Transport Tie Down

\*TRRO2 must be ordered with seat pan, seat cushion and Personal Back™.  
\*TRRO2 includes installed transport brackets, wheelchair anchored pelvic belt and locking pins.  
\*TRRO2-FD must be drop shipped to Freedom Designs for seating system installation. No Freedom Designs seating system available on-chair. Includes installed transport brackets, wheelchair anchored pelvic belt and locking pins.  
\*TRBKTS2 Includes installed Transport Brackets.

### VENT TRAY OPTION

<sup>6,7</sup>  AVENT Ventilator Tray (each).....295.00 E1029  
 ATBOX Ventilator Battery Box (each).....97.00

## ARM & BACK SELECTION

### Arm Heights

<b>AT99N</b> Range adjusts in 1" increments	<b>6.5" - 12.5"</b>
<b>AT55D</b> <b>AT55F</b> Range adjusts in 1" increments	<b>8" - 14"</b>
<b>TAD</b> <b>TAF</b>	<b>Low Range: 8" - 12"</b> <b>High Range: 10" - 16"</b>
All measurements from top of seat rail to top of armrest pad.	

### ARM TYPE

AT99N Cantilever Non Locking Bolt On Arms (pair).....No Charge  
 AT55D Cantilever Locking Desk Lgth. Bolt On Arms (pair).....195.00 E0973  
 AT55F Cantilever Locking Full Lgth. Bolt On Arms (pair).....195.00 E0973  
 TAD T-Arms, Desk Adjustable Height (pair).....225.00 E0973  
 TAF T-Arms, Full, Adjustable Height (pair) .....225.00 E0973

### T-ARM HEIGHT

L8-12 Low Range 8" - 12" .....No Charge  
 H10-16 High Range 10" - 16" .....No Charge E2618

### FOOTNOTES

1. Comes with back upholstery and adjustable angle rigidizing handle bar.
2. DPSEAT, drop mounted seat pan accommodates 2" lower seat to floor heights. Use a cushion not less than 3" thick.
3. TRRO2, TRRO2-FD, TRBKTS2 not available with Heavy Duty Reclining Back.
4. Not available with Reclining Back.
5. Transport Tie Down options not retrofittable to existing models, not field serviceable.
6. Available with 18.00 Seat to Floor Height only.
7. Available with 16" width chair only.

**Note: All specifications and dimensions are approximate.**

**\*Additional options available please contact Customer Service.**

# Invacare® Spree GT - Pediatric Growable Tilt-in-Space Wheelchair

## Accessories and Options

### CASTER & REAR SELECTIONS

#### BACK HEIGHT

- BH20 Back Height 20".....No Charge
- BH24 Back Height 24".....No Charge

#### FIXED/FOLD DOWN BACK

- FXBACK Non-Fold Down Back Canes.....**Standard**
- FLDBACK Fold-Down Back Canes.....No Charge

### FRAME FINISH SELECTIONS

#### CONTEMPORARY COLORS

- 24P Wet Black.....No Charge
- 71P Silver Vein.....No Charge
- 60P Silver Metallic.....No Charge
- 115P Black Prism.....No Charge
- 105P Midnight Blue.....No Charge

#### TOO HOT COLORS

- 30P Sunny Yellow.....No Charge
- 116P Sunkissed Orange.....No Charge
- 117P Tangerine Red.....No Charge
- 61P Electric Red.....No Charge
- 118P Fruit Fusion.....No Charge
- 119P Bubblegum Pink.....No Charge
- 120P Grape Madness.....No Charge

#### WAY COOL COLORS

- 62P Electric Blue.....No Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic Blue.....No Charge
- 123P Celestial Purple.....No Charge
- 124P Vivid Violet.....No Charge
- 125P Emerald Green.....No Charge
- 126P Shamrock Green.....No Charge
- 127P Grasshopper Green.....No Charge
- 104P Electric Teal.....No Charge

### SEAT-TO-FLOOR HEIGHT SELECTIONS (MUST WRITE IN SEAT TO FLOOR TO COMPLETE ORDER)

Seat-To-Floor.....No Charge

Front and Rear Seat-to-Floor Heights				
*Measured From Floor To Top Of Seat Rail				
Rear Wheel	4"	5"	6"	6"
	4SPLT	5SOFT 5ALUM	AT902	1259
12"	16.00"	16.00" 16.50"	16.00" 16.50" 17.00"	16.00" 16.50" 17.50" 18.00"
18"	14.50" 15.00" 15.50" 16.00"	14.50" 15.00" 15.50" 16.00"	15.50" 16.00"	16.00"
20"	15.00" 15.50" 16.00"	15.00" 15.50" 16.00" 16.50"	15.50" 16.00" 16.50" 17.00"	16.00" 16.50"
22"	16.00"	16.00" 16.50"	16.00" 16.50" 17.00"	16.00" 16.50" 17.50" 18.00"
24"	N/A	N/A	17.00"	17.50" 18.00"

All Dimensions Are Approximate

#### CASTER OPTIONS

- 4SPLT 4" Composite Urethane Light Up Casters (pair)....No Charge
- 5SOFT 5" Semi Pneumatic Casters (pair).....60.00
- 5ALUM 5" Aluminum Urethane Casters (pair).....92.00
- AT902 6" Composite Urethane 3-Spoke Casters (pair)....No Charge
- 1259 6" X 2" Semi Pneumatic Casters (pair).....60.00

#### WHEEL OPTIONS

- 12COM 12" Composite Wheels (pair).....**Standard**
- 18COM 18" Composite Wheels (pair).....No Charge
- 20COM 20" Composite Wheels (pair).....No Charge
- 22COM 22" Composite Wheels (pair).....No Charge
- 24COM7 24" Composite 7-Spoke Wheels (pair).....No Charge

#### AXLE TYPE

- A1577 Quick Release Axles (pair).....No Charge
- PAXLE Permanent Axles (pair).....**Standard**

#### TIRE OPTIONS

- PNEU Pneumatic Tires (pair).....60.00 K0067
- PNEUFF Flat Free Inserts (pair).....65.00 K0064
- URETH Urethane (Treaded) Tires (pair).....No Charge

Knee-To-Heel Range		
Front Rigging Type	Min.	Max
60°	13.50	17.50
60° MFX	9.00	12.50
70°	13.50	16.50
70° Tapered	12.75	16.25
70° MFX	5.75	11.75
90°	4.50	7.50
Elevating	14.00	20.00

Range adjust in 1" increments  
Measurement taken from top of seat rail to back of footrests.

#### HANDRIM OPTIONS

- AHDRM Aluminum Handrims (pair).....No Charge
- 1512 Plastic Coated Aluminum Handrims (pair).....72.00
- L8702 Omit Handrims.....Less 15.00

#### WHEEL LOCKS

- PULLOCK Pull To Lock Wheel Locks (pair).....No Charge
- PUSHLOCK Push To Lock Wheel Locks (pair).....No Charge
- EXTLOCK Push To Lock Wheel Locks w/6" Bolt-On Extensions (pair).65.00 E0961
- FOOTHUB Foot Operated Hublock.....200.00
- ATHUB Cane Mount Hub Lock.....50.00
- USEHUB Seat Rail Mount Hub Lock.....50.00

### ACCESSORIES & WHEEL LOCKS

#### BACK ACCESSORIES

- 1483 Non-Adjustable Rigidizing Handle Bar.....85.00
- AASH Adjustable Angle Rigidizing Handle Bar.....235.00
- ADJRGD Adjustable Width Rigidizing Bar.....85.00

#### FOOTNOTES

1. Refer to Seat-to-floor chart.
2. Not available with urethane tires.
3. Removable from back canes.
4. Non-Removable from back canes.

**Note: All specifications and dimensions are approximate.**

**\*Additional options available please contact Customer Service.**

**Invacare® Spree GT - Pediatric Growable Tilt-in-Space Wheelchair  
Accessories and Options**

**POSITIONING STRAPS**

**POSITIONING STRAPS**

E0978	<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap (each).....	50.00
E0978	<input type="checkbox"/>	1515	Push Button Buckle Seat Positioning Strap (each).....	50.00
E0978	<input type="checkbox"/>	1321BK	Hook & Loop Seat Positioning Strap (each).....	50.00
E0960	<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap (each).....	100.00
E0960	<input type="checkbox"/>	9515	Push Button Chest Positioning Strap (each).....	100.00
E0960	<input type="checkbox"/>	7321BK	Hook & Loop Chest Positioning Strap (each).....	100.00

**<sup>1</sup> HEADREST OPTIONS**

E0955	<input type="checkbox"/>	HR02S	Large Adj. Headrest System (Fixed Mount).....	284.00
E0955/E1028	<input type="checkbox"/>	CP175	Stealth Headrest 14", Multi Axis Hardware.....	305.00
E0955/E1028	<input type="checkbox"/>	CP275	Stealth Headrest 10", Multi Axis Hardware.....	305.00
E0955/E1028	<input type="checkbox"/>	CP180	Stealth Headrest 14", TWB Hardware.....	395.00
E0955/E1028	<input type="checkbox"/>	CP280	Stealth Headrest 10", TWB Hardware.....	395.00

**MISCELLANEOUS ACCESSORIES**

<input type="checkbox"/>	TLL	Tilt Lock Out Levers (pair).....	No Charge
<sup>2</sup> <input type="checkbox"/>	TLC	Tilt Lock Out Stops (pair).....	30.00
<input type="checkbox"/>	1209	Tire Pump.....	35.00
<sup>3</sup> <input type="checkbox"/>		Touch-Up Paint.....	20.00

**FRONT RIGGING OPTIONS**

**FRONT RIGGING**

E0990	<input type="checkbox"/>	70LIFTNT	70° Lift Off Non-Taper Footrests (pair).....	95.00
	<input type="checkbox"/>	70LIFTMFX	70° Lift Off MFX Footrests (pair).....	95.00
K0053	<input type="checkbox"/>	90LIFTMFX	90° Lift Off MFX Footrests (pair).....	95.00
	<input type="checkbox"/>	SWING1319	60° Swingaway Footrests (pair).....	No Charge
	<input type="checkbox"/>	SWING914	60° Swingaway Footrests MFX (pair).....	No Charge
	<input type="checkbox"/>	70NTAPER	70° Swingaway Footrests (pair).....	No Charge
E0990	<input type="checkbox"/>	70MFX	70° Swingaway Footrests MFX (pair).....	No Charge
E0990	<input type="checkbox"/>	70STAPER	70° Tapered Footrests (pair).....	No Charge
	<input type="checkbox"/>	90SWING	90° Swingaway Footrests (pair).....	No Charge
	<input type="checkbox"/>	AT5044	Elevating Legrests (pair).....	255.00
K0053	<input type="checkbox"/>	BILAT	Bi-Lateral Contracture Footrests (pair).....	120.00
K0108	<input type="checkbox"/>	ATL4	Center Pivot Articulating Legrests (pair).....	320.00
	<sup>4</sup> <input type="checkbox"/>	110CP	Contracture Platform.....	280.00

**FOOTPLATE TYPE**

<sup>5</sup> <input type="checkbox"/>	1651	Flip-Up Composite Footplates (pair).....	No Charge	
	A1562	Aluminum Footplates (pair).....	59.00	
K0041	<sup>6</sup> <input type="checkbox"/>	A1564	Extra Large Aluminum Footplates (pair).....	59.00
K0040	<sup>7</sup> <input type="checkbox"/>	AT5543	Adjustable Angle Flip-Up Footplates (pair).....	165.00

**ADJUSTABLE ANGLE FOOTPLATES**

(FOR BI-LATERAL CONTRACTUAL FOOTRESTS ONLY)

K0040	<input type="checkbox"/>	AAF325R	Adj. Angle Footplate Right 3.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF425R	Adj. Angle Footplate Right 4.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF525R	Adj. Angle Footplate Right 5.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF625R	Adj. Angle Footplate Right 6.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF725R	Adj. Angle Footplate Right 7.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF325L	Adj. Angle Footplate Left 3.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF425L	Adj. Angle Footplate Left 4.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF525L	Adj. Angle Footplate Left 5.25" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF625L	Adj. Angle Footplate Left 6.26" W (each).....	83.00
K0040	<input type="checkbox"/>	AAF725L	Adj. Angle Footplate Left 7.25" W (each).....	83.00

**FOOTREST/LEGREST ACCESSORIES**

K0038	<input type="checkbox"/>	1337	Calf Strap (each).....	24.00
K0039	<input type="checkbox"/>	9342	"H" Strap (each).....	56.00
	<input type="checkbox"/>	IMPACTG	Impact Guards (pair).....	51.00
	<sup>8</sup> <input type="checkbox"/>	ALPT	Longer Pivot And Slide Tubes (pair).....	84.00
	<input type="checkbox"/>	ALPT3	3" Extension For 90° Footrests (pair).....	31.00

**HEEL LOOP/TOE LOOP**

<input type="checkbox"/>	A1600BK	Heel Loops With Ankle Straps (pair).....	51.00	E0951
<input type="checkbox"/>	HLOOP	Heel Loops Without Ankle Straps (pair).....	51.00	E0951
<input type="checkbox"/>	1345	Toe Loops (pair).....	39.00	E0952

**FOOTNOTES**

1. Headrest Options only available on-chair w/Infinity and Personal Backs.
2. Additional pair of tilt out stops, one pair standard with chair.
3. Touch-up paints are not available in all colors.
4. Not available on seat widths less than 16".
5. Not available on 90SWING, 70MFX, 90LIFTMFX or 70LIFTMFX front riggings.
6. Not available with 70STAPER.
7. Dorsi/plantar angle adjustable.
8. ALPT is available on 70° non-taper and tapered footrests on 16"- 22" wide chairs, using a composite or adjustable angle footplates. ALPT is not available on 70MFX or extra large footplates.

**Note: All specifications and dimensions are approximate.  
\*Additional options available please contact Customer Service.**

## Invacare® Spree GT - Pediatric Growable Tilt-in-Space Wheelchair Accessories and Options

### Specifications

- |                         |  |
|-------------------------|--|
| ♦ Frame Style           | Tilt-in-Space                          |
| ♦ Seat Widths           | 12" - 18"                              |
| ♦ Seat Depths           | 12" - 19"                              |
| ♦ Back Height           | 20" and 24" (See Back Height Category) |
| ♦ Adjustable Angle Back | 0°, 10°, 20°, 30°                      |
| ♦ Camber Angles         | 0°                                     |
| ♦ Tubing O.D.           | 1" Reg<br>Seat rail 1"                 |
| ♦ Weight Limit          | 200 lbs                                |
| ♦ Weight of Chair       | 26 lbs*<br>38 lbs**                    |
| ♦ Warranty              | Lifetime on frames and crossbraces     |

\*Base chair weight without accessories

\*\*16"x16", 90SWING, AT5543, AT902, 12COM, PNEU, Anti-tippers,  
No Seat Pan, No Arms, No Wheel Locks, Subtract 4.9 lbs to get weight  
w/o 90SWING and AT5543.

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING**

**INFINITY® SEAT AND BACK CUSHIONS  
SEE PRICE LIST FORM # 01-351.**

**CONVENTIONAL PRODUCTS  
SEE PRICE LIST FORM # 95-82.**

**ABSOLUTE™, LAP TRAYS AND  
PERSONAL BACK PRODUCTS  
SEE PRICE LIST FORM # 01-207.**

### Basic Chair Includes:

- ♦ Spree GT Tilt-in-Space Frame
- ♦ Transport Tie Down Options
- ♦ Seat Width (W12 - W16)
- ♦ Seat Depth (D12 - D16)
- ♦ Arm Type (AT99N)
- ♦ Back Height Selection
- ♦ Fixed or Fold Down Back
- ♦ Frame Finish
- ♦ Caster Options (4SPLT or AT902)
- ♦ Wheel Options
- ♦ Axle Type Selection
- ♦ URETH Urethane (treaded) Tires
- ♦ Handrim Options (AHDRM)
- ♦ Front Rigging Selection (70NTAPER, SWING1319, SWING914,  
90SWING, 70STAPER, or 70MFX)
- ♦ 1651 Flip Up Composite Footplates
- ♦ TLL Tilt Out Levers (1 pair)
- ♦ TLC Tilt Lock Out Stops (1 pair)
- ♦ Wheel Locks (PUSHLOCK or PULLLOCK)
- ♦ Anti-tippers

### Available Through Our Specials Department

Adjustable Clamp-on Sector Blocks (pair)	SM-031007103339
--	-----------------

Invacare Corporation      [www.invacare.com](http://www.invacare.com)



**Yes, you can.®**

**USA**  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

Invacare, Yes, you can., Infinity and the Medallion Design are registered trademarks  
Absolute is a trademark of Invacare Corporation. Specifications and prices are  
subject to change without notification.

©2006 Invacare Corporation

Form No. 04-143 Rev 07/06  
P2004\_143I\_SPREE GT