

Invacare® Spree XT

Price List and Order Form Effective September 15, 2003 Revised October 18, 2005



For ease of ordering, contact Customer Service Toll Free at:
 800-333-6900 Fax: 800-678-4682
 To order Literature: 800-828-6282

www.invacare.com

Yes, you can.™

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME TYPE

E1234 SPRXT Spree XT Folding Tilt2,895.00

TRANSPORT TIE DOWN

- TRRO Transport Ready OptionNo Charge
- TRBKTS Transport Brackets & Pelvic Belt.....No Charge
- TROMIT Less Transport Ready Option.....Less 150.00

Transport Tie Down

- * TRRO must be ordered with Invacare KSS system and includes transport brackets and wheelchair anchored pelvic belt
- * TRBKTS includes installed brackets and wheelchair anchored pelvic belt
- * TRRO/TRBKTS not available with 2" seat rail extension, side guards, AT5044 Elevating legrests

SEAT WIDTH

- W10 10" Wide (Growable to 12").....No Charge
- W11 11" Wide (Growable to 12").....No Charge
- W12 12" Wide (Growable to 14").....No Charge
- W13 13" Wide (Growable to 14").....No Charge
- W14 14" Wide (Growable to 16").....No Charge
- W15 15" Wide (Growable to 16").....No Charge
- W16 16" Wide (Growable to 18").....No Charge
- W17 17" Wide (Growable to 18").....No Charge
- W18 18" WideNo Charge

SEAT DEPTH RANGE

- FT1214 Short Seat Rail (12" - 14" Deep).....No Charge
- FT1218 Long Seat Rail (12" - 18" Deep).....No Charge

SEAT RAIL EXTENSION

- SRE2 2" Seat Rail Extension Installed25.00

SEAT DEPTH SETTING

- D12 Seat Depth 12"No Charge
- D13 Seat Depth 13"No Charge
- D14 Seat Depth 14"No Charge
- D15 Seat Depth 15"No Charge
- D16 Seat Depth 16"No Charge
- D17 Seat Depth 17"No Charge
- D18 Seat Depth 18"No Charge
- D19 Seat Depth 19"No Charge
- E2203 D20 Seat Depth 20"300.00

Arm Heights

AT99N	6" - 9"
Range adjusts in 1" increments	
AT55D	9" - 12"
AT55F	
Range adjusts in 1" increments	
TAD	Low Range: 8" - 12"
TAF	High Range: 10" - 16"
All measurements from top of seat rail to top of armrest pad.	

ARM TYPE

- AT99N Cantilever Non Locking Bolt On Arms (pair).....No Charge
- AT55D Cantilever Locking Desk Lgth. Bolt On Arms (pair).....195.00 E0973
- AT55F Cantilever Locking Full Lgth. Bolt On Arms (pair).....195.00 E0973
- TAD T-Arms, Desk Adjustable Height225.00 E0973
- TAF T-Arms, Full, Adjustable Height.....225.00 E0973
- L8401 Omit Arms.....Less 20.00
- L8405 Omit Arms Pads.....Less 10.00

T-ARM HEIGHT

- L8-12 Low RangeNo Charge
- H10-16 High Range.....No Charge

ARM ATTACHMENT

- BOLT Bolt-On.....No Charge
- CLAMP Clamp-On.....40.00

FIXED BACK HEIGHT

- BH20 20" Height.....No Charge
- BH24 24" Height.....No Charge

FOOTNOTES

1. A seating system must be purchased for use with this wheelchair.
2. Transport Tie Down options not retrofittable to existing models, not field serviceable.
3. No seat pan option is offered with this wheelchair.
4. Setting requires selection of FT1218 OR FT1214 and SRE2.
5. Setting requires selection of FT1218 and SRE2.
6. Not available with AT99N.
7. Back does not come with upholstery.

NOTE: All specifications and dimensions are approximate.

Invacare® Spree XT

Accessories and Options

FRAME FINISH

CONTEMPORARY COLORS

- 24P Wet Black.....No Charge
- 71P Silver Vein.....No Charge
- 60P Silver Metallic.....No Charge
- 115P Black Prism.....No Charge
- 105P Midnight Blue.....No Charge

TOO HOT COLORS

- 30P Sunny Yellow.....No Charge
- 116P Sunkissed Orange.....No Charge
- 117P Tangerine Red.....No Charge
- 61P Electric Red.....No Charge
- 118P Fruit Fusion.....No Charge
- 119P Bubbegum Pink.....No Charge
- 120P Grape Madness.....No Charge

WAY COOL COLORS

- 62P Electric Blue.....No Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic Blue.....No Charge
- 123P Celestial Purple.....No Charge
- 124P Vivid Violet.....No Charge
- 125P Emerald Green.....No Charge
- 126P Shamrock Green.....No Charge
- 127P Grasshopper Green.....No Charge
- 104P Electric Teal.....No Charge

CASTER OPTIONS

- AT900 3" Roller Casters (pair).....38.00
- 4SOFT 4" Semi Pneumatic Casters (pair).....60.00
- 4SPLT 4" Composite Urethane Light Up Casters (pair).....No Charge
- 5SOFT 5" Semi Pneumatic Casters (pair).....60.00
- AT901 5" Composite Urethane Casters (pair).....64.00
- AT902 6" Composite Urethane 3 Spoke Casters (pair).....No Charge
- 6PNEU 6" Composite Pneumatic 6 Spoke Casters (pair).....64.00
- 1259 6" X 2" Semi Pneumatic Casters (pair).....60.00
- AT903 8" X 1 1/4" Composite Urethane Casters (pair).....78.00
- AT1292FF Flat Free Inserts for 8" X 1 1/4" Composite Pneu Casters (pair).....60.00
- 1292 8" X 1 1/4" Composite Pneumatic Casters (pair).....60.00

SEAT-TO-FLOOR HEIGHT

(MUST WRITE IN SEAT TO FLOOR TO COMPLETE ORDER)

_____ Seat-To-Floor.....No Charge

Seat-To-Floor Heights							
*Measured From Floor To Top of Seat Rail							
Rear Wheel	3"	4"	5"	6"	6"	8"	Tilt Range
	AT900	4SPLT 4SOFT	AT901 5SOFT	AT902 6PNEU	1259	AT903	
12"	16.00 16.50 17.00	15.50 16.00 16.50 17.00	16.00 16.50 17.00 17.50	17.00 17.50	17.50	N/A	45°
20"	N/A	N/A	17.50	17.50 18.00 18.50	17.50 19.00	18.50	45°
22"	16.50 17.00	16.50 17.00	16.50 17.00 17.50 18.00	17.00 17.50	17.50	N/A	30°
24"	N/A	N/A	17.50	17.50 18.00 18.50	17.50 19.00	18.50	30°

All Dimensions Are Approximate

WHEEL OPTIONS

- 12COM 12" Composite Wheels (pair).....No Charge
- 20COM 20" Composite Wheels (pair).....No Charge
- 22COM 22" Composite Wheels (pair).....No Charge
- 24COM7 24" Composite 7 Spoke Wheels (pair).....No Charge
- 20SP 20" Spoke Wheels (pair).....No Charge
- 22SP 22" Spoke Wheels (pair).....No Charge
- 24SP 24" Spoke Wheels (pair).....No Charge

TIRE OPTIONS

- PNEU Pneumatic Tires (pair).....60.00 K0067
- PNEUFF Flat Free Inserts (pair).....65.00 K0064
- URETH Urethane (Treaded) Tires (pair).....No Charge
- PRIMO Primo Tires (pair).....60.00 K0067
- KIKB KIK® Black Tires (pair).....60.00

AXLE TYPE

- A1577 Quick Release Axles (pair).....No Charge
- PAXLE Permanent Axles (pair).....No Charge

HANDRIM OPTIONS

- ALUMPC Aluminum Anodized Handrims (pair).....No Charge
- A1512K Black Plastic Coated Handrims (pair).....72.00 K0059
- AT1224-8 Plastic Coated w/8 Oblique Projections (pair).....245.00 E0967
- AT1225-8 Plastic Coated w/8 Vertical Projections (pair).....245.00 E0967
- AT1225-12 Plastic Coated w/12 Vertical Projections (pair).....245.00 E0967
- L8702 Omit Handrims.....Less 15.00

MOUNTING HANDRIMS (SPOKE ONLY)

- IN Mounting Inside Handrim.....No Charge
- OUT Mounting Outside Handrim.....No Charge

Knee-To-Heel Range		
Front Rigging Type	Min.	Max
60°	13.50	17.50
60° MFX	9.00	12.50
70°	13.50	16.50
70° Tapered	12.75	16.25
70° MFX	5.75	11.75
90°	4.50	7.50
Elevating	14.00	20.00

Range adjust in 1" increments
Measurement taken from top of seat rail to back of footrests.

FOOTNOTES

1. Not available on 12COM, 20SP, 22SP and 24SP.
2. Only available on 24SP.
3. Only available on 24COM7 and 24SP.
4. Not available with 12COM.
5. Only available on 24COM7.
6. Decreases overall width by 3/4".

NOTE: All specifications and dimensions are approximate.

Invacare® Spree XT

Accessories and Options

FRONT RIGGING

<input type="checkbox"/>	SWING1319	60° Swingaway Footrests (pair)No Charge
<input type="checkbox"/>	SWING914	60° Swingaway Footrests MFX (pair)No Charge
<input type="checkbox"/>	70NTAPER	70° Swingaway Footrests (pair)No Charge
¹ <input type="checkbox"/>	70STAPER	70° Tapered Footrests (pair)No Charge
<input type="checkbox"/>	70MFX	70° Swingaway Footrests MFX (pair)No Charge
<input type="checkbox"/>	90SWING	90° Swingaway Footrests (pair)No Charge
E0990 ² <input type="checkbox"/>	AT5044	Elevating Legrests (pair)255.00

FOOTPLATE TYPE

K0040 ³ <input type="checkbox"/>	AT5543	Adjustable Angle Flip-Up Footplates (pair)165.00
<input type="checkbox"/>	1651	Flip-Up Composite Footplates (pair)No Charge
⁴ <input type="checkbox"/>	A1562	Aluminum Footplates (pair)59.00

HEEL LOOP/TOE LOOP

E0951 <input type="checkbox"/>	A1600BK	Heel Loops With Ankle Straps (pair)51.00
E0951 <input type="checkbox"/>	HLOOP	Heel Loops Without Ankle Straps (pair)51.00
E0952 <input type="checkbox"/>	1345	Toe Loops (pair)39.00

FOOTREST/LEGREST ACCESSORIES

K0038 <input type="checkbox"/>	1337	Calf Strap (each)24.00
⁵ <input type="checkbox"/>	IMPACTG	Impact Guards (pair)51.00
⁶ <input type="checkbox"/>	ALPT	Longer Pivot And Slide Tubes (pair)No Charge
<input type="checkbox"/>	ALPT2	Foot Support Extension 2" for 90° (pair)25.00
<input type="checkbox"/>	ALPT3	Foot Support Extension 3" for 90° (pair)31.00
K0039 <input type="checkbox"/>	A9342	"H" Strap (each)51.00

POSITIONING STRAPS

E0978 <input type="checkbox"/>	1515	Push Button Buckle Seat Positioning Strap (each)50.00
E0978 <input type="checkbox"/>	1321BK	Hook & Loop Seat Positioning Strap (each)40.00
E0960 <input type="checkbox"/>	9515	Push Button Chest Positioning Strap (each)100.00
E0960 <input type="checkbox"/>	7321BK	Hook & Loop Chest Positioning Strap (each)100.00

BACK ACCESSORIES

<input type="checkbox"/>	BP02	Back Pack50.00
<input type="checkbox"/>	1483	Non-Adjustable Rigidizing Handle Bar Standard
<input type="checkbox"/>	AASH	Adjustable Angle Rigidizing Handle Bar235.00
<input type="checkbox"/>	RIGB	Flip Down Rigidizer Bar85.00
<input type="checkbox"/>	SHOMIT	Omit Non Adjust Rigidizer Handle BarLess 85.00

SEAT ACCESSORIES

<input type="checkbox"/>	SP-01	Seat Pouch35.00
--------------------------	-------	------------	------------

7 ARM ACCESSORIES

<input type="checkbox"/>	1671-SM	Rigid Side Guard, Small (pair)122.00
<input type="checkbox"/>	1671-LG	Rigid Side Guard, Large (pair)122.00
<input type="checkbox"/>	1671-SMH	Quad Rigid Side Guard Small (pair)122.00
<input type="checkbox"/>	1671-LGH	Quad Rigid Side Guard Large (pair)122.00

WHEEL LOCKS

⁸ <input type="checkbox"/>	9691	Push to Lock Clamp On Wheel Locks (pair)No Charge
<input type="checkbox"/>	A1567	Pull To Lock Clamp On Wheel Locks (pair)No Charge
E0974 ⁸ <input type="checkbox"/>	A1598	Hill Holder Clamp On Wheel Locks (pair)77.00
E0961 <input type="checkbox"/>	A1576	Wheel Lock Extensions (8") (pair)39.00
⁹ <input type="checkbox"/>	ATHUB	Cane Mount Hublock (attendant)50.00
⁹ <input type="checkbox"/>	USEHUB	Seat Rail Mount Hublock50.00

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	3421	3/4" To 1" Clamp Assembly (each)21.00
K0065 ¹⁰ <input type="checkbox"/>	A1307	Spoke Protectors (For 20",22",24",26") (pair)90.00
¹¹ <input type="checkbox"/>	TLC	Tilt Lock Out Collars (pair)35.00
<input type="checkbox"/>		Touch Up Paint20.00

FOOTNOTES

1. 16" - 18" widths only.
2. 14" - 18" widths only.
3. Dorsi/plantar angle adjustable.
4. Not available on 90SWING and 70MFX.
5. Not available with 90SWING, 70MFX, SWING914 and AT5044.
6. 2" longer slide tube used on composite and adjustable angle footplates only.
Only available with 60°, 70°, and 70° Tapered.
7. Not available with TAD and TAF.
8. Not available with 22" and 24".
9. Only available with 12COM & 20COM rear wheels
10. Not available with projection handrims.
11. Additional pair of tilt lock out collars, one pair standard with chair.

NOTE: All specifications and dimensions are approximate.

Invacare® Spree XT Accessories and Options

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING**

**INFINITY® SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 01-351.**

**CONVENTIONAL PRODUCTS
SEE PRICE LIST FORM # 95-82.**

**ESSENTIAL, LAP TRAYS, AND
PERSONAL BACK™ PRODUCTS
SEE PRICE LIST FORM # 01-207.**

Weight Limitation 150 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Specifications

♦ Frame Style	Folding Tilt
♦ Seat Widths	10" - 18"
♦ Seat Depths	12" - 20" (with seat rail extension)
♦ Back Angle Adj.	-4° to 16°
♦ Tubing Size	1"
♦ Tilt Range	0° - 45°
♦ Overall Chair Width	Seat width plus 7.50" (12, 20) Seat width plus 8.75" (22, 24)
♦ Warranty	Lifetime on frames and crossbraces

Basic Chair Includes:

- ♦ Spree XT Folding Tilt
- ♦ TRRO or TRBKTS Transport Tie Down Option
- ♦ Seat Width Selection
- ♦ Seat Depth Setting (D12 - D19)
- ♦ Arm Type (AT99N)
- ♦ Bolt On Arm Attachment
- ♦ Fixed Back Height Selection
- ♦ Frame Finish
- ♦ Caster Options (2SPLT or AT902)
- ♦ Seat-to-Floor Height Selection
- ♦ Wheel Options
- ♦ URETH - Urethane (Treaded) Tires
- ♦ Axle Type (A1577 or PAXLE)
- ♦ ALUMPC - Aluminum Anodized Handrims
- ♦ Front Rigging Selection - SWING1319, SWING914, 70NTAPER, 70STAPER, 70MFX or 90SWING
- ♦ I651 - Flip-Up Composite Footplates (pair)
- ♦ ALPT - Longer Pivot and Side Tubes
- ♦ I483 - Non-Adjustable Rigidizing Handle Bar
- ♦ Wheel Lock Selection (9691 or A1567)
- ♦ Anti-Tippers (pair)
- ♦ Tilt Trigger Lock-Out Levers
- ♦ Tilt Lock Out Collars (1 pair)

Specials Department

For your customization requests, please contact Customer Service for details and prices.

Growth Kit Available

Invacare Corporation www.invacare.com



Yes, you can.™

USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

Invacare, Infinity and the medallion design are registered trademarks and Yes, you can. and Personal Back are trademarks of Invacare Corporation. KIK is a registered trademark of KIK Technology, Inc.
Specifications and prices are subject to change without notification.
©2005 Invacare Corporation

Form No. 03-016 Rev. 10/05
P2003_016G_Spree XT