

Invacare Matrix PSP Cushion



The Invacare Matrix PSP (Posture Seat Polymer) cushion features an anatomically contoured shape with a unique waffled ischial relief recess combined with a polymer overlay for extra skin protection and shear reduction.

Our Matrix PSP Cushions are available in Standard adult sizes (16W (41cm), 17W (43cm), 18W (46cm), 20W (51cm)) and Pediatric sizes (14W (36cm)).

Weight Capacity: Standard = 300 lb. (136kg); Pediatric = 150 lb. (68kg)

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit

Configured model numbers are: **PSP**wwdd



i Refer to the **Invacare Matrix PS Cushion User Manual** (1171735) for important Safety Information and detailed instructions on proper cushion Usage/Installation and Care & Maintenance.



TRD0503 Rev A (supplement to 1171735)

Invacare Matrix PSP Cushion



The Invacare Matrix PSP (Posture Seat Polymer) cushion features an anatomically contoured shape with a unique waffled ischial relief recess combined with a polymer overlay for extra skin protection and shear reduction.

Our Matrix PSP Cushions are available in Standard adult sizes (16W (41cm), 17W (43cm), 18W (46cm), 20W (51cm)) and Pediatric sizes (14W (36cm)).

Weight Capacity: Standard = 300 lb. (136kg); Pediatric = 150 lb. (68kg)

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit

Configured model numbers are: **PSP**wwdd



i Refer to the **Invacare Matrix PS Cushion User Manual** (1171735) for important Safety Information and detailed instructions on proper cushion Usage/Installation and Care & Maintenance.



TRD0503 Rev A (supplement to 1171735)