

Pegasus Pro

Feel the Freedom



Comfort, safety and reliability
with a powerful drive performance



Yes, you can.®

Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From ergonomic steering controls to the adaptable cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

Quality



Market-leading build quality with exceptional attention to detail.

Safety



Compliant with the latest European Medical Device Standards.

Comfort



Adjustable design to provide lasting comfort.



Feel in control with our **Pegasus Pro**

Introducing our Pegasus Pro, a high performing scooter designed for your comfort, safety and reliability wherever you go.

Ideal for longer journeys, the Pegasus Pro features advanced suspension and large wheels to ensure you drive in comfort and safety even across uneven ground,

around corners or over obstacles. Incredibly simple to steer and start, the Pegasus Pro puts you in total control.





Improved driving control

- Ergonomic steering for added control
- Designed for all-terrain performance
- Anti-splash guards



Safety in mind

- Automatic speed reduction for safe cornering
- Two-step disengaging lever to prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Designed for **improved traction around corners** and safe, all-terrain performance





Built around comfort

- High quality cushioned seat
- Adjustable steering column
- Adjustable seat height
- Flip-up armrests and swivel seat makes getting on and off easy



Adjustable steering for greater comfort



Swivel seat for easy transfer



High quality lighting



Two-step disengaging lever to prevent free-wheeling



Add to your scooter

Tailor the Orion® Pro to your taste and needs. Turn over to explore the optional extras available.

Add to your scooter

We have a wide range of optional extras available, so you can customise your Pegasus Pro scooter. Ask your retailer to discover the full range.



Rear basket



Side & rear view mirrors



Front & rear bumper



Mobile phone holder



Crutch holder



Seatbelt



Cupholder



Interchangeable shroud panels

Technical data

For more information about this product, including the user manual, please visit our local website.



Backrest width ▼



Seat depth ▼



Seat to floor height ▼



Backrest height ▼



Driving unit width ▼



Total length ▼

3 wheels	510 mm	460 mm	Min. 440 mm Max. 510 mm	630 mm	650 mm	1300 mm
4 wheels	510 mm	460 mm	Min. 440 mm Max. 510 mm	630 mm	650 mm	1320 mm



Empty weight ▼



The heaviest part ▼



Max. User weight ▼



Battery capacity ▼



Turning diameter ▼



Max. Climbable obstacle height ▼

3 wheels	126 kg	62 kg	127 kg	2 x 12 V/75 Ah (C20) AGM	2500 mm	100 mm
4 wheels	136 kg	62 kg	160 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm



Ground clearance ▼



Max. Safe slope ▼



Motor capacity ▼



Max. Driving distance ▼



Max. Driving speed ▼

3 wheels	100 mm	10°	10 km/h: 240W/600W ¹ 12.8 km/h: 550W/1600W ¹ 15 km/h: 550W/1600W ¹	52 km	10 / 12.8 / 15 km/h
4 wheels	100 mm	10°	6 km/h: 240W/500W ¹ 10 km/h: 240W/600W ¹ 12 km/h: 550W/1500W ¹ 12.8 km/h: 550W/1600W ¹ 15 km/h: 550W/1600W ¹	52 km	6 / 10 / 12 / 12.8 / 15 km/h

¹ Peak performance

Colours*



Marble White



Alium Blue



Diamond Silver



Garnet Red



Jade Green



Luxor Beige



Crystal Copper



Carmine Red



Ionic Chrome



Electric Blue



Canary Yellow



Black Vinyl

* Scooter colours may vary slightly. Please check with your retailer.

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2026



Yes, you can.®